•• PRODUCT DESCRIPTION

Dental paste for permanent filling of root canals of permanent and temporary teeth based on zinc oxide and eugenol. In the humid environment of the root canal a zinc eugenolate – insoluble salt is forming after mixture.

•• COMPOSITION

- Main components: zinc oxide;
- Additional components: hydrocortisone acetate, dexamethasone phosphate, thymol iodide, barium sulfate, bismuth oxynitrate, magnesium stearate.

•• THE MAIN PROPERTIES AND CHARACTERISTICS

- Does not dissolve after hardening (in case of root-canal obturation with gutta-percha pins).
- Does not shrink over time.
- Has increased radiopacity.
- Can be easily inserted into the root canal and, if necessary, may be easily removed from the root canal.
- Is tolerated by periapical tissues.
- Dexamethasone phosphate and hydrocortisone acetate provide postoperative sensitization decrease and anti-inflammatory effect.
- Thymol iodide causes the antimicrobial effect of the paste.
- The antiseptic effect of the paste continues for several hours after being introduced into the canal in order to sterilize organic residues in the canal that may remain after the pulp extirpation.
- The antiseptic and anti-inflammatory effect of the paste stops after the hardening.

•• INDICATIONS FOR USAGE

Permanent root canal filling of permanent and temporary teeth and sealing of root canals in combination with gutta-percha. Paste is compatible with all obturation techniques:

- single pin technique;
- lateral condensation;
- vertical condensation and compaction;
- thermal obturation.

•• MODE OF APPLICATION

After the finishing of the medicinal and instrumental treatment of the root canals of the tooth, determination of the working length of the canals, in the absence of contraindications for filling (obturation) of the root canal, it is necessary to isolate the working field (using cotton rolls or rubber dam). Before filling, it is recommended to treat root canals with an antiseptic solution and dry them with paper points or other hygroscopic products.

Mixing method: powder «Dex-&-Dent» is mixed with eugenol in a volume ratio of 5:1, resulting in a paste with the consistency of thick sour cream.

The paste may be injected into the root canals using canal fillers (Lentulo) or reamers (files, spreaders) according to the technique of filling (obturation) of the tooth root canals chosen by the doctor. After the end of the procedure, it is recommended to perform a control X-ray examination and, if necessary, improve the canal filling.

- **Working time: minimum 6 hours at room temperature.**
- **Hardening in the root canal: from 1 till 12 hours.**

•• CONTRAINDICATIONS

- Presence in the patient's history of allergic reactions to any component of the material or their intolerance
- Improper use of the product.

•• INTERACTIONS WITH OTHER DENTAL MATERIALS

Zinc oxide-eugenol root canal filling paste is not compatible with dental composite filling materials. The components of the paste can interfere with the polymerization of composites.

•• PRECAUTIONARY MEASURES

- Do not swallow! Causes burns.
- Avoid getting the paste and its components on the surface of mucous membranes and skin.
- The use of protective gloves and face and eye protection is recommended.
- In case of complications or patient feeling unwell, it is necessary to immediately consult a doctor, who, if possible, should be shown this instruction.

•• SIDE-EFFECTS

An allergic reaction to the material is possible if it or its components come into contact with mucous membranes or skin.

•• SPECIAL WARNINGS

The material is used only after complete instrumental and medical treatment of the root canals of the teeth.

It is not recommended to allow the material to exit the apex of the tooth root (the root canal is filled strictly to its working length)

Paste removed from the cortex canal can be resorbed for a long time and causes irritation of the periapical tissues.

Like all powder forms, the material is electrostatic and tends to form granules during transport or storage. Therefore, before use, the bottle should be shaken several times.

In case of contact with eyes, oral mucosa or skin, rinse with copious amounts of water and consult a doctor.

•• WARRANTY

The manufacturer warrants the quality of manufactured products. The adverse events inflicted by violation of user manual, storage conditions and other events inflicted by non-stipulated usage of the material are not the subjects of warranty.

The customer is responsible for determination of suitability of this product for user's application. Warranty conditions: the product does not comply with requirements declared by manufacturer. In this case the manufacturer replaces the defective material within warranty period.

•• LIMITATION OF LIABILITY

The manufacturer's liability is limited by only cases stipulated by direct legislation of the country.

•• STORAGE AND SHELF LIFE

- Store the material in a dry place at a temperature of 4 - 25 °C.
- Do not expose the product to direct sunlight.
- Inadequate storage conditions will reduce the terms of use and may lead to a deterioration of the material properties.
- Warranty storage period is 2 years.

•• RECYCLING

Dispose of the medical device in accordance with local / regional / national / international legal requirements.

•• PHARMACEUTICAL FORM

- Jar of 10 g powder
- Jar of 20 g powder

The preparation can be supplied with Jen-Eugenol liquid bottle 6 or 10 ml for paste mixing.

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